GEARFIL SEDIO 2 Dental Universal Self-Etch Adhesive

Flow Chart Sheet







Kuraray Noritake Dental Inc.

16<mark>21 Sakaz</mark>u, Kurashiki, Okayama 710-0801, Japan 001 1561R987R

INDICATIONS

CASE

[1] Direct restorations using light-cured composite resin

Direct restorations

- [2] Cavity sealing as a pretreatment for indirect restorations*
- [3] Treatment of exposed root surfaces*
- [4] Treatment of hypersensitive teeth*
- [5] Intraoral repairs of fractured restorations

Intraoral repairs

- [6] Post cementation using a dual- or self-cured composite resin*
- [7] Core build-ups using a light-, dual- or self-cured core material*
- [8] Cementing inlays, onlays, crowns, bridges and veneers using a composite resin cement



Cementing

*Please refer to the IFU for [2], [3], [4], [6] and [7] of Indications.

Table: Dental curing unit and curing time


Туре	Light source	Light Intensity	Light-curing time
Halogen	Halogen lamp	More than 400 mW/cm ²	10 seconds
LED	Blue LED*	800-1400 mW/cm ²	10 seconds
		More than 1500 mW/cm ²	5 seconds

The effective wavelength range of each dental curing unit must be 400-515 nm.

^{*}Peak of emission spectrum: 450-480nm

inlays, onlays, crowns, bridges and veneers using a composite resin cement

Clean and dry the tooth surface, and then trial fit the prosthetic restoration



CLEARFIL SE BOND 2 Kuraray Noritake Dental Inc.



*Before applying PRIMER, selective enamel etch by phosphoric acid is an option.







*Light-curing the mixture is an option for optimum performance.



of fractured restorations















*Refer to the table for light-curing time.



*Use an opaque resin (e.g. CLEARFIL ST OPAQUER) to mask metal color.

using light-cured composite resin

Follow the standard procedures for isolation, moisture control, cavity preparation and pulp protection



*Before applying PRIMER, selective enamel etch by phosphoric acid is an option.









*Refer to the table for light-curing time.

