

product OVERVIEW



Implantati TiUltra

NobelParallel Conical Connection TiUltra

Platforma	Implantat Ø	Dužina*						
		7 mm	8.5 mm	10 mm	11.5 mm	13 mm	15 mm	18 mm
	3.75 mm	300295	300296	300297	300298	300299	300300	300301
	4.3 mm	300302	300303	300304	300305	300306	300307	300308
	5.0 mm	300309	300310	300311	300312	300313	300314	300315
	5.5 mm	300316	300317	300318	300319	300320	300321	



* Važno: stvarna dužina Implantata je 0.5 mm kraća od napisane dužine.

NobelActive TiUltra

Platforma	Implantat Ø	Dužina*						
		7 mm	8.5 mm	10 mm	11.5 mm	13 mm	15 mm	18 mm
	3.0 mm			300240	300241	300242	300243	
	3.5 mm		300244	300245	300246	300247	300248	300249
	4.3 mm		300250	300251	300252	300253	300254	300255
	5.0 mm		300256	300257	300258	300259	300260	300261
	5.5 mm	300262	300263	300264	300265	300266	300267	



* Važno: stvarna dužina Implantata je 0.5 mm kraća od napisane dužine.

NobelReplace Conical Connection TiUltra

Platforma	Implantat Ø	Dužina*				
		8 mm	10 mm	11.5 mm	13 mm	16 mm
	3.5 mm	300337	300338	300339	300340	300341
	4.3 mm	300342	300343	300344	300345	300346
	5.0 mm	300347	300348	300349	300350	300351



* Važno: stvarna dužina Implantata je 0.6 mm duža od napisane dužine.

Pokrovni vijak

3.0	NP	RP	WP
36775	36649	36650	37812



NobelActive PureSet

Svrdla, screw taps i bone mills se naručuju odvojeno



NobelParallel CC PureSet

Svrdla, screw taps i bone mills se naručuju odvojeno



NobelReplace CC PureSet

Bone mills se naručuju odvojeno

All-in-one Set: NobelActive i NobelParallel CC kirurški instrumenti

Zatražite ponudu za All-in-one Set

Implantati TiUnite

NobelParallel Conical Connection TiUnite*

Platforma	Implantat Ø	Dužina						
		7 mm	8.5 mm	10 mm	11.5 mm	13 mm	15 mm	18 mm
	3.75 mm	37963	37964	37965	37966	37967	37968	37969
	4.3 mm	37970	37971	37972	37973	37974	37975	37976
	5.0 mm	37977	37978	37979	37980	37981	37982	37983
	5.5 mm	37984	37985	37986	37987	37988	37989	



NobelActive TiUnite

Platforma	Implantat Ø	Dužina						
		7 mm	8.5 mm	10 mm	11.5 mm	13 mm	15 mm	18 mm
	3.0 mm			36769	36770	36771	36772	
	3.5 mm		35221	34125	34126	34127	34128	35215
	4.3 mm		35223	34131	34132	34133	34134	35219
	5.0 mm		35225	34137	34138	34139	34140	35220
	5.5 mm	37806	37807	37808	37809	37810	37811	



NobelReplace Conical Connection TiUnite

Platforma	Implantat Ø	Dužina				
		8 mm	10 mm	11.5 mm	13 mm	16 mm
	3.5 mm	36699	36700	36701	36702	36703
	4.3 mm	36704	36705	36707	36708	36709
	5.0 mm	36710	36711	36712	36713	36714



TiUnite
na vratu

NobelReplace Conical Connection PMC TiUnite* (djelomično strojno obrađen vrat)

Platforma	Implantat Ø	Dužina				
		8 mm	10 mm	11.5 mm	13 mm	16 mm
	3.5 mm	37284	37285	37287	37288	37289
	4.3 mm	37290	37291	37292	37293	37294
	5.0 mm	37295	37296	37297	37298	37299



0,75 mm
strojno
obrađenog
vrata

Pokrovni vijak

3.0	NP	RP	WP
36775	36649	36650	37812



Oznake implantata

Protetske oznake
Protetske komponente
kodirane su bojom za točno
i brzo prepoznavanje.



Platform
3.0



Narrow
Platform
NP



Regular
Platform
RP



Wide
Platform
WP



DTX Studio Implantat

Ključ za uspješno liječenje

Trodimenzionalno dijagnosticirajte i planirajte Implantate na temelju anatomije pacijenata i protetskih zahtjeva pomoću programa DTX Studio Implantat. Za složenije slučajeve mogućnost postavljanja Implantata uz potpuno vođenje Nobel Guided Surgery



Doživotna garancija Nobel Biocare

Na sve Nobel Biocare Implantate uključujući prefabricirane protetske komponente.

Za više informacija posjetite nobelbiocare.com/warranty

Kontaktirajte NOVODENT i tražite posebne ponude za DTX Studio Implantat www.novodent.hr

* Uključen pokrovni vijak.

Kirurški instrumenti

Svrkla i nareznice

Precision Drill	
	36118

Guide Drill	
	35426

NobelActive & NobelParallel Conical Connection

Twist Drills	7-10 mm	7-15 mm	10-18 mm
Ø 1.5 mm		31278	
Ø 2.0 mm	32296	32297	32299



Twist Step Drills	7-10 mm	7-15 mm	10-18 mm
Ø 2.4/2.8 mm	32260	32261	32262
Ø 2.8/3.2 mm	37873	34638	34639
Ø 3.2/3.6 mm	32263	32264	32265
Ø 3.8/4.2 mm	32275	32276	32277
Ø 4.2/4.6 mm	37874	34582	34583
Ø 4.2/5.0 mm	37875	37876	



Screw Taps NobelParallel CC	7-10 mm	7-13 mm	7-15 mm	7-18 mm
Ø 3.75 mm		37990		37991
Ø 4.3 mm		37992		37993
Ø 5.0 mm		37994		37995
Ø 5.5 mm	37996		37997	



Cortical Drills NobelParallel CC	
Ø 3.75 mm	38000
Ø 4.3 mm	38001
Ø 5.0 mm	38002
Ø 5.5 mm	38003



Screw Taps NobelActive	
Ø 3.0 mm	36816
Ø 3.5 mm	36236
Ø 4.3 mm	36237
Ø 5.0 mm	36238
Ø 5.5 7.0-10.0 mm	37871
Ø 5.5 11.5-15.0 mm	37872



NobelReplace Conical Connection

Drill with Tip Tapered	
Ø 2.0 mm	36117



Drills	8 mm	10 mm	11.5 mm	13 mm	16 mm
NP Ø 3.5 mm	32075	29367	36113	29368	29369
RP Ø 4.3 mm	32076	29370	36114	29371	29372
WP Ø 5.0 mm	32077	29373	36115	29374	29375



Dense Bone Drills	13 mm	16 mm
NP Ø 3.5 mm	29377	29378
RP Ø 4.3 mm	29380	29381
WP Ø 5.0 mm	29383	29384



Screw Taps	
NP Ø 3.5 mm	36717
RP Ø 4.3 mm	32090
WP Ø 5.0 mm	32091




Kirurške komponente za sve CC Implantate

Implantat Drivers	28 mm	37 mm
3.0	36773	36774
NP	36718	36719
RP	36720	36721
WP	37859	37860



Bone Mills and Guides	Bone Mill with Guide	Bone Mill Guide
3.0 Ø 4.0 mm	37861	37862
NP Ø 4.4 mm	37863	37865
NP Ø 5.2 mm	37864	37865
RP Ø 5.2 mm	37866	37868
RP Ø 6.2 mm	37867	37868
WP Ø 6.7 mm	37869	37870



32068 All-on-4® Guide
Kirurška šablona za All-on-4 tehniku

Nadogradnje za cijeljenje

	Visina	3.0			NP			RP			WP	
		3.0mm	5.0mm	7.0mm	3.0mm	5.0mm	7.0mm	3.0mm	5.0mm	7.0mm	3.0mm	5.0mm
	Healing Abutment											
	Ø 3.2mm	36794	36795	36796								
	Ø 3.6mm				36639	36640	36867	36643	36644	36872		
	Ø 3.8mm	36797	36798	36799								
	Ø 5.0mm				36641	36642	36868	36645	36646	36873	37813	37814
	Ø 6.0mm							36647	36648	36874		
	Healing Abutment Bridge¹											
	Ø 4.0mm				36864	36865	36866					
	Ø 5.0mm							36869	36870	36871		
	Ø 6.0mm										37817	37818

¹Napomena

Koristiti za Universal Base Non-Engaging i ostale Non-Engaging protetske komponente za radove na vijak na razini implantata
Potrebno koristiti otisne kapice Bridge Open Tray i ScanBody za nadomjestak više zubi (bez oznake SA)

	Healing Abutment Anatomical PEEK*										6.0×7.0mm	7.0×8.0mm
											37819	37820

	Slim Healing Abutment		37669	37670		37666	37665		37667	37668		
---	------------------------------	--	-------	-------	--	-------	-------	--	-------	-------	--	--

Implantat Drivers za Slim Healing i Temporary Abutments

	Implantat Driver CC 3.0										37713	
	Implantat Driver CC NP										37677	
	Implantat Driver CC RP										37678	

Kirurške šablone - lokalno printanje

Free Surgical Template (DTX Studio Implant)

10419

Pilot Drill procedura

Guided Pilot Drill Sleeve	1.5	2.0	
	300438	300440	
Mounting Tool Pin Guide	1.5	2.0	1.5 / 2.0
	300442	300443	300444
Fully guided Sleeve	NP	RP	WP / 6.0
	32754	32765	32766

Fully Guided procedura

Fully guided Sleeve	NP	RP	WP / 6.0
	32754	32765	32766
Guided Anchor Pin	1.5	1.5 short sharf	
	30900	34761	
Guided Cylinder Tool CC Implant	NP	RP	WP
Ø 3.5mm	37172		
Ø 4.3mm		37173	
Ø 5.0mm		37174	
Ø 5.5mm			37950



NobelActive Guided PureSet

Svrkla, screw taps i bone mills se naručuju odvojeno



NobelParallel CC Guided PureSet

Svrkla, screw taps i bone mills se naručuju odvojeno



NobelReplace CC Guided PureSet

Bone mills se naručuju odvojeno

Protetska rješenja

Finalne nadogradnje

1 zub		Više zubi		Potpuna bezubost	
Vijak	Cementiranje	Vijak	Cementiranje	Vijak	Mobilno
Universal Base (engaging/1-zub)	Esthetic Abutment				
On1 Universal Abutment (engaging/1-zub)	On1 Esthetic Abutment (engaging/1-zub)	On1 Universal Abutment (non-engaging/most)	On1 Esthetic Abutment (engaging/most)	Universal Base MUA Straight Non-Engaging (non-engaging/most)	Locator R-Tx®
		Universal Base Straight Non-Engaging (non-engaging/most)			
		Universal Base MUA Straight Non-Engaging (non-engaging/most)			

Privremene nadogradnje

Temporary Snap Abutment Engaging (Single-unit)*		NP	RP	WP	
	Visina 1.5 mm	38760	38761	38762	
	Visina 3.0 mm	38847	38848	38849	
Temporary Abutment Engaging (Single-unit)*		3.0	NP	RP	WP
	Visina 1.5 mm	36779	36663	36664	37823
	Visina 3.0 mm				37824
Temporary Abutment Non-Engaging (Bridge)*			NP	RP	WP
	Visina 1.5 mm		36661	36662	37825
	Visina 3.0 mm				37826
Slim Temporary Abutment		3.0	NP	RP	
	Visina 6.5 mm	37675	37671	37673	
	Visina 7.5 mm	37676	37672	37674	
Temporary Abutment Anatomical PEEK*					WP
	Visina 6.0×7.0 mm				37821
	Visina 7.0×8.0 mm				37822


Otisne kapice Razina implantata

Closed Tray	3.0	NP	RP	WP
	Visina 13.0mm	13.0mm	9.0mm	13.0mm
	Ø 3.3 mm	36801		
	Ø 3.6 mm	36538	36541	36540
	Ø 3.8 mm	36803		
	Ø 5.0 mm	36539	36543	36542
	Ø 6.0 mm		36545	36544
	Ø 6.5 mm			37853
				37852
Open Tray	3.0	NP	RP	WP
	Visina 14.0mm	10.0mm	14.0mm	10.0mm
	Ø 3.2 mm	36800		
	Ø 3.6 mm	36258	36260	36263
	Ø 3.8 mm	36802		
	Ø 5.0 mm	36259	36261	36265
	Ø 6.0 mm		36267	36266
	Ø 6.5 mm			37857
				37856
Bridge Open Tray ¹		NP	RP	WP
	Visina	12.0 mm	12.0 mm	12.0 mm
		36930	36931	37858

¹ Napomena

Koristiti za Universal Base Non-Engaging i ostale Non-Engaging protetske komponente za radove na vijak na razini implantata
Potrebno koristiti Healing Abutment Bridge i ScanBody za nadomjestak više zubi (bez oznake SA)

Konzistentno kodiranje bojom
Protetske komponente
kodirane su bojom za točno
i brzo prepoznavanje.



Esthetic Abutment

			3.0			NP			RP			WP		
	Esthetic Abutment*	Visina	1.5 mm	3.0 mm	4.5 mm	1.5 mm	3.0 mm	4.5 mm	1.5 mm	3.0 mm	4.5 mm	6.0×7.0 mm	7.0×8.0 mm	
			36782	36783	36814	36665	36666	36249	36669	36671	36251	37827	37828	

	Esthetic Abutment 15° *		36784	36785	36815	36667	36668	36250	36672	36673	36252			
--	-------------------------	--	-------	-------	-------	-------	-------	-------	-------	-------	-------	--	--	--

			NP			RP			WP		
	Procera® Titan Abutment*	Visina	1.5 mm	3.0 mm	3.0 mm kut 20°	1.5 mm	3.0 mm	3.0 mm kut 20°	3.0 mm kut 15°		
		Ø 4,5mm	✓	✓	✓	✓	✓	✓			
		Ø 5,5mm	✓	✓	✓	✓	✓	✓			
		Ø 6,5mm				✓	✓		✓		

Universal Base

			3.0	NP		RP		WP	
	Universal Base*	Visina		1.5 mm	3.0 mm	1.5 mm	3.0 mm	1.5 mm	3.0 mm
				38213	38216	38214	38217	38215	38218

	Universal Base Straight Non-Engaging*1	Visina		1.5 mm		1.5 mm		1.5 mm	
				301362		301363		301364	

!Napomena za Non-Engaging!

Potreban OmniGrip zavijač

Potrebno koristiti Healing Abutment Bridge, otisne kapice Bridge Open Tray
i ScanBody za nadomjestak više zubi (bez oznake SA)

Laboratorijski implantati - Lab Replica

Laboratorijski Implantati (na razini implantata)	3.0	NP	RP	WP
 Lab Replica - gipsani model	36791	36697	36698	37879

	Lab IOS Replica - printani model Single-unit & Multiple-unit restorations**		PMA-NBA35-1	PMA-NBA43-1	PMA-NBA60-1
--	---	--	-------------	-------------	-------------

** Za ugradnju koristiti instrumente: PMA-AIP-2 / PMA-AIS-5 / PMA-AIT-1

Instrumenti za ugradnju Lab Replica u printani model		
	Insertion Pin	PMA-AIP-2
	Insertion Screw	PMA-AIS-5
	Pliers	PMA-AIT-1

Scan Body za intraoralno i desktop skeniranje (ordinacija i laboratorij)

IOS Scan Body (Clinic & Lab)		3.0	NP	RP	WP	Kit
	Single-unit restorations Nadomjestak 1 zuba	IO 2B-D SA	IO 2B-A SA	IO 2B-B SA	IO 2B-C SA	IO 2B SA Kit
	Multiple-unit restorations Nadomjestak više zubi		IO 2B-A	IO 2B-B	IO 2B-C	IO 2B Kit
	Multi-unit Abutment		IO 2C-A	IO 2C-A	IO 2C-B	IO 2C Kit

Instrumenti za postavljanje IOS Scan Body		
	Dugi	C13048
	Kratki	C13485

Multi-unit Abutments Xeal & Universal Base MUA

	NP	RP	WP
	Multi-unit Abutment Xeal		
	Visina 1.5 mm	300171	300172 300173
	Visina 2.5 mm	300174	300175 300176
	Visina 3.5 mm	300177	300178 300179
	Visina 4.5 mm		300180

	17° Multi-unit Abutment Xeal*		
	Visina 2.5 mm	300181	300182 300183
	Visina 3.5 mm	300184	300185 300186

	30° Multi-unit Abutment Xeal*		
	Visina 3.5 mm	300187	300188
	Visina 4.5 mm	300189	300190

Universal Base Multi-unit Abutment Straight* / **			
	NP / RP	Visina 0.3 mm	301365
	WP	Visina 0.3 mm	301366

****potreban OmniGrip Mini zavijač; protetski vijak MUA Omni Grip Mini**

Multi-unit Aligning Instrument	
	300161

Healing Cap Multi-unit (titan)	(2/pkg)	Ø 5.0	Ø 6.0	Wide
	Visina 4.1 mm	300162	300164	300166
	Visina 5.5 mm	300163	300165	300167

Healing Cap Multi-unit (plastic)	
	(5/pkg) 29064
	(wide 1/pkg) 31146

Temporary Snap Coping Multi-unit*	
	38915

Temporary Coping Multi-unit*	
	29046

Otisne kapice	
	Impression Coping Open Tray Multi-unit 29089
	Impression Coping Closed Tray Multi-unit 38924

IOS Scan Body & Laboratorijski Implantati (na razini Multi Unita MUA)	3.0	NP	RP	WP
	Lab Replica MUA - gipsani model		38918 (1/ pkg)	38918 (1/ pkg) 31162
			38919 (5/pkg)	38919 (5/pkg)
	Lab IOS Replica MUA - printani model**		PMA-MUA45-1	PMA-MUA45-1 PMA-MUA60-1
	IOS Multi-unit Abutment		IO 2C-A	IO 2C-A IO 2C-B

**** Za ugradnju koristiti instrumente: PMA-AIP-2 / PMA-AIS-5 / PMA-AIT-1**

Abutment Retrieval Rescue Set

Instrumenti za rješavanje komplikacija sa nadogradnjama

29161	Handle for Machine Instruments	37502	Abut Screw Retr Reverse Drill RP/WP/6.0
37480	Handle for Impl Rescue Collar & Guides	37503	Abutment Screw Retr Instrument 3.0/NP
37485	Rescue Drill Guide CC 3.0	37504	Abutment Screw Retr Instrument RP/WP/6.0
37486	Rescue Drill Guide CC NP	37505	Abutment Screw Remover 3.0
37487	Rescue Drill Guide CC RP	37506	Abutment Screw Remover NP
37497	Screw Tap Repair M1.4	37507	Abutment Screw Remover RP/WP/6.0
37498	Screw Tap Repair M1.6	37517	Screw Tap Repair M2.0
37499	Screw Tap Repair M1.8	37933	Rescue Drill Guide CC WP
37500	Screw Tap Repair M2.5	301344	Abutment Screw Retrieval Set Organizer
37501	Abutment Screw Retr Reverse Drill 3.0/NP		

Implant Retrieval Rescue Set

Instrumenti za rješavanje komplikacija sa implantatima

37473	Impl Retr Instr CC 3.0 31mm
37474	Impl Retr Instr CC NP & Ext Hex WP 22mm
37475	Impl Retr Instr CC RP & Tri-Ch WP 22mm
37477	Trephine Drill 3.2/4.0mm
37480	Handle for Impl Rescue Collar & Guides
37927	Impl Retr Instr CC WP 22mm
37928	Trephine Drill 3.8/4.6mm
37929	Trephine Drill 4.4/5.2mm
37930	Trephine Drill 5.2/6.2mm
37931	Trephine Drill 5.6/6.6mm
37932	Trephine Drill 6.2/7.0mm
301343	Implant Retrieval Set Organizer

Lokatori

Locator Abutments

Locator® Abutment			
	1.0 mm	REF02069	REF02075
	2.0 mm	REF02070	REF02076
	3.0 mm	REF02071	REF02077
	4.0 mm	REF02072	REF02078
	5.0 mm	REF02073	REF02079
	6.0 mm	REF02074	REF02080

Locator® coping	4/pkg	20/pkg
	REF08505	

Locator Process kit	
	REF08519-10

Locator® Abutment Replica	4/pkg
	Ø 4 mm REF08530
	Ø 5 mm REF08516

Instrumenti za Locator Abutments		
	Latch Tip Torque Wrench Driver	REF08914
	Core Tool	REF08393

Locator R-Tx® Abutments u kompletu sa Process kit

Locator R-Tx® Abutment	NP	RP	WP
	Visina 1.0mm	REF30506-01	REF30507-01 REF30508-01
	Visina 2.0mm	REF30506-02	REF30507-02 REF30508-02
	Visina 3.0mm	REF30506-03	REF30507-03 REF30508-03
	Visina 4.0mm	REF30506-04	REF30507-04 REF30508-04
	Visina 5.0mm	REF30506-05	REF30507-05 REF30508-05
	Visina 6.0mm	REF30506-06	REF30507-06

Locator R-Tx® Impression Coping	4/pkg
	REF30017-01

Processing components Locator R-Tx®	(4/pkg)
	Processing insert REF30012-01
	Denture attachment processing assembly REF30013-01
	Processing spacer REF30018-01
	Zero retention insert REF30001-01
	Low retention insert REF30002-01
	Medium retention insert REF30003-01
	High retention insert REF30004-01

Locator R-Tx® Abutment replicas	4/pkg	20/pkg
	Ø 3.35 mm	REF30014-01
	Ø 4.0 mm	REF30015-01
	Ø 5.0 mm	REF30016-01

Instrumenti za Locator R-Tx®		
	Hex Latch Driver Torque Wrench 0.050" (1.25 mm)	REF08008
	Retention Insert Tool	REF30021-01

Nobel Procera®

Make the right choice
for your patients
NobelProcera®
Restorative Solutions



NobelProcera nadogradnje i krunice na implantatima

- ✓ Odabir individualne nadogradnje od titana, cirkona bez Ti adaptera i cirkonske krunice u punoj formi sa anguliranim pristupom.
- ✓ 100% cement-free protetska rješenja na vijak.
- ✓ Velika protetska fleksibilnost korištenjem angulacijskog kanala za vijak.
- ✓ Brz i praktičan rad sa Omnigrip instrumentima i vijcima

NobelProcera strukture mostova i prečke

- ✓ Titanske strukture na implantatima i na MultiUnit abutmentima za 2-14 članova
- ✓ Cirkonske strukture na implantatima i na MultiUnit abutmentima za 2-14 članova
- ✓ Protetske strukture za mostove na vijak
- ✓ Prečke u različitim oblicima i formama

Individualne nadogradnje i krunice

Titanska



nadogradnja po implan-
tatu

Cirkonska



bez Ti adaptera



sa Ti adapterom

nadogradnja sa
korekcijom kuta
0-25° ASC



sa Ti adapterom

krunica u punoj anatomske
formi sa korekcijom kuta
0-25° FCZ

Individualni mostovi na implantatima i MUA

NobelProcera Implant Bridge Ti*



titanska
struktura

Implant level
(max. 14 članova)

MUA level
(max. 14 članova)

*bez mogućnosti korekcije kuta

NobelProcera Implant Bridge



cirkonska
struktura

most u punoj formi
sa mogućnosti
korekcije kuta 0-25°

Individualne prečke na implantatima

NobelProcera Bar Ti



struktura prečke
od titana

2-3 implantata

4-5 implantata

NobelProcera® protetsko rješenje za sve vaše slučajeve

Od nedostatka jednog zuba do potpunih bezubosti, kod slučajeva na zubima ili implantatima, NobelProcera ima inovativno i znanstveno dokazano protetsko rješenje za sve potrebe pacijenata.



Doktor



Slanje otiska
standardno
ili digitalno

Dizajn sa
DTX Studio™
Lab

NobelProcera
Scan & Design
Service

Dizajn sa
3Shape®



NobelProcera
proizvodnja



Dobivate
individualno rješenje

Protetski instrumenti

PureSet protetski set

PureSet Basic protetski set			301233
		UniGrip Screwdriver Machine 20 mm	29151
		UniGrip Screwdriver Machine 30 mm	29153
		Multi-unit Screwdriver Machine 21 mm	29158
		OmniGrip Screwdriver Machine 20 mm*	37379
		OmniGrip Screwdriver Machine 30 mm*	37381
		Manual Torque Wrench Prosthetic	29165

Klinički i laboratorijski vijci i zavijač

Esthetic Abutment, Universal Base, Temporary Abutment, NobelProcera
Individualne titanske nadogradnje, strukture mosta od titana (na implantatu i MUA) i prečke od titana

Zavijač: UniGrip

		3.0	NP	RP/WP
	Klinički vijak	37890	37891	37892
	Laboratorijski vijak	36805	37894	37895 (5/pkg)

UniGrip™ Screwdriver Machine		
	20 mm	29151
	25 mm	29152
	30 mm	29153
	35 mm	29154

Multi-unit Abutment i protetske radove na MUA / NobelProcera
Strukture mosta od titana i prečke na MUA

Zavijač: UniGrip

		NP	RP/WP
	Vijak Multi-unit kutni Abutment	36892	37893
	Protetski vijak Multi-unit Abutment	29285	29285
	Laboratorijski protetski vijak Multi-unit Abutment	29287 (5/pkg) 38420 (1/pkg)	29287 (5/pkg) 38420 (1/pkg)

UniGrip™ Screwdriver Manual		
	20 mm	29148
	28 mm	29149
	36 mm	29150

Universal Base Non engaging Straight, NobelProcera

Individualne cirkonske nadogradnje sa korekcijom kuta 0-25° ASC (Angulated Screw Channel) i Individualne krunice u anatomske forme sa korekcijom kuta 0-25° FCZ
Strukture mosta od cirkona sa korekcijom kuta 0-25° na razini implantata

Zavijač: OmniGrip

		NP	RP/WP
	OmniGrip klinički vijak	37367	37606
	OmniGrip laboratorijski vijak	37374	37607

OmniGrip™ Screwdriver Machine		
	20 mm	37379
	25 mm	37380
	30 mm	37381
	35 mm	37382

OmniGrip™ Screwdriver Manual		
	20 mm	37376
	28 mm	37377
	36 mm	37378

Universal Base Multi-unit Abutment Straight NobelProcera

Strukture mosta od cirkona sa korekcijom kuta 0-25° na MUA

Zavijač: OmniGrip Mini

		NP/RP	WP
	Protetski vijak MUA OmniGrip Mini	301203	301200
	Laboratorijski protetski vijak MUA OmniGrip Mini	301201	301202


OmniGrip™ Mini Screwdriver Machine		
	20 mm	300852
	28 mm	300853
	36 mm	300854

OmniGrip™ Mini Screwdriver Manual		
	20 mm	300855
	36 mm	300856

On1

On1 Base Xeal incl. handle		NP	RP	WP
	1.75 mm	300191	300193	300195
	2.5 mm	300192	300194	300196

Pakiranje sa On1 kliničkim vijkom



Klasični postupak

On1 Healing Cap		NP	RP	WP
	1.5 mm	38695	38697	38699
	2.5 mm	38696	38698	38700

Digitalni postupak

On1 IOS Healing Cap incl. handle		NP	RP	WP
	4 mm			38739
	4.5 mm	38735	38737	
	5 mm			38738
	6 mm	38734	38736	

Pakiranje sa On1 protetskim vijkom

- IOS Position Locator za DTX Studio, 3Shape i Exocad
- Desktop Position Locator za 3Shape i Exocad

On1 Clinical Screw		NP	RP	WP
	1.75 mm	38744	38746	38748
	2.5 mm	38745	38747	38749

On1 Base Screwdriver		
	Machine 21 mm	38752
	Manual 25 mm	38753

On1 Temporary Abutment incl. handle		NP	RP	WP
	Engaging (Single-unit)	38701	38703	38705
	Non-engaging (Bridge)	38702	38704	38706

Pakiranje sa On1 protetskim vijkom

On1 Universal Abutment		NP	RP	WP
	Engaging (Single-unit)			
	0.3 mm	38707	38709	38711
	1.25 mm	38708	38710	38712
	Non-engaging (Bridge)			
	0.3 mm	300361	300363	300365
	1.25 mm	300362	300364	300366

Pakiranje sa On1 protetskim vijkom i On1 plastičnim cilindrom

On1 Base Replica (za gipsani model)	NP	RP	WP
	38731	38732	38733

On1 IOS Replica (za printani model)	NP	RP	WP
	38740	38741	38742

On1 Laboratory Prosthetic Screw	NP	RP	WP
	38751	38751	38751

On1 IOS Replica Insertion Tool	NP	RP	WP
	38201	38201	38201

On1 Impression Coping		NP	RP	WP
	Closed Tray	38725	38726	38727
	Open Tray	38728	38729	38730
	Open Tray	300367	300368	300369
	Non-engaging (Bridge)			

On1 Esthetic Abutment Titanium		NP	RP	WP
	0.3 mm	38713	38714	38715

Pakiranje sa On1 protetskim vijkom

On1 Prosthetic Screw	NP	RP	WP
	38750	38750	38750

Instrumenti za imedijatnu temporizaciju

Apical Drill	
	38853

Implantat protection analog/drill guide	NP	RP	WP
	38855 (5/pkg)	38856 (5/pkg)	38857 (5/pkg)

	Drill Guide Multi-unit	38917
---	------------------------	-------

	Protection Analog Multi-unit	29123 (5/pkg)
---	------------------------------	------------------

	Držac za Protective Analog	29122
--	----------------------------	-------

NobelPearl™



Protetska nadogradnja

Ravna (dvije visine vrata) i kutna dostupna za implantate 3.5 mm (NP), 4.2 mm (RP) i 5.5 mm (WP).

Vicarbo vijak

Posebno ojačani vijak sa karbonskim vlaknima.

Digitalni workflow

NobelPearl Ceramic Base: Individualno oblikovana nadogradnja za protetske radove na vijak. Prilagođeno za rad u DTX, exocad i 3Shape softveru.

NobelPearl Tapered (sa pokrovnim vijkom)		NP	RP	WP
	8 mm	300752	300755	300759
	10 mm	300753	300756	300760
	12 mm	300754	300757	300761
	14 mm		300758	

NobelPearl Healing Abutments		NP	RP	WP
	Ø 4mm 3 mm	300702		
	4 mm	300705		
	Ø 5mm 3 mm		300708	
	4 mm		300709	
	Ø 6mm 3 mm			300710
	4 mm			300712

NobelPearl Abutments (sa vijkom)		NP	RP	WP
	Straight 1 mm	300663	300664	300665
	2 mm	300666	300667	300668
	15° 1 mm	300669	300670	300671

NobelPearl Ceramic Base (sa vijkom)		NP	RP	WP
	Engaging	301263	301264	301265
	Non-Engaging	301285	301286	301287

NobelPearl Temporary Abutment ¹		NP	RP	WP
		300713	300714	300715

¹ Uključen Temporary Clinical Screw

NobelPearl Docklocs Abutments		2 mm	3 mm	4 mm
		301252	301253	301254

NobelPearl Impression Coping		NP	RP	WP
	Open Tray	300716	300717	300718
	Closed Tray	300719	300720	300721
	Closed Tray Long	300722	300723	300724

NobelPearl Position Locator ³		NP	RP	WP
		301246	301247	301248

³ Uključen Position Locator Screw

NobelPearl Implant Replica ⁴		NP	RP	WP
		301249	301250	301251

⁴ Bez Digital Implant Replica Placer

NobelPearl Tapered Drill Stop		NP	RP	WP
		300766	300767	300768

300772 NobelPearl Prosthetic Kit Inter-X

301307 Abutment Screw Retrieval Kit

301267 NobelPearl Tapered Surgery Kit

GalvoSurge®

A NEW APPROACH
to
implant cleaning

Aparat za čišćenje
površine implantata
Inovativno rješenje kod
periimplantitisa
2 minute period čišćenja
Jednostavno za korištenje

GalvoSurge
Control Unit GS 1000
GVS1002
- dental implant cleaning system



GalvoSurge Dental Implant
Cleaning Set
GVS1008



Torque kontrola

Implantati


NobelActive® TiUnite 3.0
NobelActive® TiUltra 3.0



Max 45 Ncm




NobelActive® TiUnite NP, RP, WP
NobelActive® TiUltra NP, RP, WP



Ostali Nobel Biocare implantati

NobelParallel® Conical Connection TiUnite
NobelReplace® ConicalConnection TiUnite
NobelReplace® ConicalConnection PMC TiUnite

NobelParallel® Conical Connection TiUltra
NobelReplace® ConicalConnection TiUltra



Protetske komponente

Healing Abutment
Slim Healing
(incl. Healing Abutment Anatomical PEEK)
Healing Cap Multi Unit
Healing CapMulti-unit Titan
Otisne kapice
On1 Healing Cap
On1 Impression Coping



Ručno

Slim Temporary Abutments
Protetski vijak
Titanium Temporary Coping
Titanium Temporary Snap Coping



17° Multi-unit Abutment Xeal
30° Multi-unit Abutment Xeal
Universal Base Multi Unit




Abutments for
NobelActive® 3.0



Max 15 Ncm



Sve ostale privremene
i finalne nadogradnje
(uključujući Temporary Abutment Anatomical PEEK)
Locator
Locator R-tx
Universal Base ENG.
Universal Base NON-ENG
Universal Base On1
On1 Base
Multi-Unit Abutment Xeal



Kompatibilni softveri



3shape

exocad

dental wings



Preuzmite biblioteke za
kompatibilne softvere
i lokalnu proizvodnju

www.nobelbiocare.com/local-production



novodent

Uvoznik & distributer

Novodent d.o.o.
info@novodent.hr
prodaja@novodent.hr
www.novodent.hr

MATULJI
Cesta dalmatinskih brigada 7
51211 Matulji • 051 701 645

ZAGREB
Zavrtnica 3
10000 Zagreb • 01 457 7846

SPLIT
Ruđera Boškovića 5B (Fregata II. kat)
21000 Split • 091 270 4712 / 091 170 1645