



Universal Base

Product overview

nobelbiocare.com/universal-base



Universal Base – Engaging

Single unit restoration



Universal Base engaging

Connection	CC			External hex			Tri-channel		
Platform	NP	RP	WP	NP	RP	WP	NP	RP	WP
Article nr. collar 1.5 mm	38213	38214	38215	38227	38228	38229	38171	38172	38173
Article nr. collar 3.0 mm	38216	38217	38218	38230	38231	38232	38174	38225	38226
Post height	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm
Diameter	Ø 4.775 mm	Ø 4.775 mm	Ø 6.515 mm	Ø 4.775 mm	Ø 4.775 mm	Ø 6.515 mm	Ø 4.775 mm	Ø 4.775 mm	Ø 6.515 mm
Screw	37891	37892	37892	29282	29283	29284	36818	29475	29475
Screwdriver	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip
Lab Screw	37894	37895*	37895*	37894	37895*	31169	31170	29293	29293
Elos Accurate™ Position Locator	IO 2B-A SA	IO 2B-B SA	IO 2B-C SA	IO 6A-A	IO 6A-B	IO 6A-C	IO 2A-A	IO 2A-B	IO 2A-C
Elos Accurate™ Position Locator Kit	IO 2B SA Kit			IO 6A Kit			IO 2A Kit		
Replica for IOS - Elos Accurate™ Analog	PMA-NBA35-1	PMA-NBA43-1	PMA-NBA60-1	PMA-BRA35-1	PMA-BRA41-1	PMA-BRA51-1	PMA-NBR35-1	PMA-NBR43-1	PMA-NBR50-1
Desktop scanbody DTX	300457	300458	300459	300460	300461	300462	300452	300453	300454
Plaster model implant analog	36697	36698	37879	31158	31159	31160	29498	29500	29502

* 5/pkg

Universal Base Engaging: Burn-out coping and clinical screw included.

Universal Base – Non-Engaging

Multiple unit restoration – Conical design

35
Uni



Universal Base non-engaging – Conical design


Connection	CC			External hex			Tri-channel		
Platform	NP	RP	WP	NP	RP	WP	NP	RP	WP
Article nr.	301235	301236	301237	301240	301241	301242	301243	301244	301245
Abutment collar height	1.8 mm	1.0 mm	1.0 mm	1.5 mm	1.0 mm	0.8 mm	1.0 mm	0.6 mm	0.6 mm
Post height	4.1 mm	4.1 mm	4.1 mm	4.1 mm	4.1 mm	4.1 mm	4.1 mm	4.1 mm	4.1 mm
Diameter	Ø 5.0 mm	Ø 5.0 mm	Ø 5.7 mm	Ø 5.0 mm	Ø 5.0 mm	Ø 6.0 mm	Ø 5.0 mm	Ø 5.0 mm	Ø 5.7 mm
Screw	37891	37892	37892	29282	29283	29284	36818	29475	29475
Screwdriver	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip
Lab Screw	37894	37895*	37895*	37894	37895*	31169	31170	29293	29293
Elos Accurate™ Position Locator	IO 2B-A	IO 2B-B	IO 2B-C	IO 6A-A	IO 6A-B	IO 6A-C	IO 2A-A	IO 2A-B	IO 2A-C
Elos Accurate™ Position Locator Kit	IO 2B Kit			IO 6A Kit			IO 2A Kit		
Replica for IOS - Elos Accurate™ Analog	PMA-NBA35-1	PMA-NBA43-1	PMA-NBA60-1	PMA-BRA35-1	PMA-BRA41-1	PMA-BRA51-1	PMA-NBR35-1	PMA-NBR43-1	PMA-NBR50-1
Desktop scanbody DTX	300470	300471	300472	300463	300464	300465	300466	300467	300468
Plaster model implant analog	36697	36698	37879	31158	31159	31160	29498	29500	29502

* 5/pkg

Universal Base non-engaging: clinical screw included.

Universal Base – Non-Engaging

Multiple unit restoration – Straight design



Universal Base non-engaging – Straight design



Connection	CC			External hex			Tri-channel		
Platform	NP	RP	WP	NP	RP	WP	NP	RP	WP
Article nr.	301362	301363	301364	301367	301368	301369	301370	301371	301372
Abutment collar height	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm
Post height	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm
Diameter	Ø 3.6 mm	Ø 4.2 mm	Ø 5.1 mm	Ø 3.9 mm	Ø 4.4 mm	Ø 5.4 mm	Ø 3.7 mm	Ø 4.5 mm	Ø 5.1 mm
Screw	37367	37606	37606	29282	29283	29284	36818	29475	29475
Screwdriver	OmniGrip	OmniGrip	OmniGrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip
Lab Screw	37374	37607	37607	37894	37895*	31169	31170	29293	29293
Elos Accurate™ Position Locator	IO 2B-A	IO 2B-B	IO 2B-C	IO 6A-A	IO 6A-B	IO 6A-C	IO 2A-A	IO 2A-B	IO 2A-C
Elos Accurate™ Position Locator Kit	IO 2B Kit			IO 6A Kit			IO 2A Kit		
Replica for IOS - Elos Accurate™ Analog	PMA-NBA35-1	PMA-NBA43-1	PMA-NBA60-1	PMA-BRA35-1	PMA-BRA41-1	PMA-BRA51-1	PMA-NBR35-1	PMA-NBR43-1	PMA-NBR50-1
Desktop scanbody DTX	300470	300471	300472	300463	300464	300465	300466	300467	300468
Plaster model implant analog	36697	36698	37879	31158	31159	31160	29498	29500	29502

* 5/pkg

Universal Base non-engaging: clinical screw included.

Universal Base – Multi unit

Multiple unit and full arch restoration – Conical design

15
Uni



Universal Base Multi unit – Conical design		
Connection	Multi-Unit	
Platform	NP / RP ¹	WP ²
Article nr.	301238	301239
Abutment collar height	0.9 mm	1.2 mm
Post height	4.1 mm	4.1 mm
Diameter	Ø 5.0 mm	Ø 5.7 mm
Screw	29285	29286
Screwdriver	Unigrip	Unigrip
Lab Screw	29287*	29287*
Elos Accurate™ Position Locator	IO 2C-A	IO 2C-B
Elos Accurate™ Position Locator Kit	IO 2C Kit	
Replica for IOS - Elos Accurate™ Analog	PMA-MUA45-1	PMA-MUA60-1
Desktop scanbody DTX	300473	300474
Plaster model implant analog	38918	31162

* 5/pkg

¹ Compatible with all Multi-unit Abutments EXCEPT Multi-unit Abutment WP ex-hex (Brånemark)

² Compatible ONLY with Multi-unit Abutments WP ex-hex (Brånemark)

Universal Base non-engaging: clinical screw included.

Universal Base – Multi unit

Multiple unit and full arch restoration – Straight design

15
O-Mini



Universal Base Multi unit – Straight design		
Connection	Multi-Unit	
Platform	NP / RP ¹	WP ²
Article nr.	301365	301366
Abutment collar height	0.3 mm	0.3 mm
Post height	4.5 mm	4.5 mm
Diameter	Ø 5.0 mm	Ø 5.9 mm
Screw	301203	301200
Screwdriver	OmniGrip Mini	OmniGrip Mini
Lab Screw	301201	301201
Elos Accurate™ Position Locator	IO 2C-A	IO 2C-B
Elos Accurate™ Position Locator Kit	IO 2C Kit	
Replica for IOS - Elos Accurate™ Analog	PMA-MUA45-1	PMA-MUA60-1
Desktop scanbody DTX	300473	300474
Plaster model implant analog	38918	31162

¹ Compatible with all Multi-unit Abutments EXCEPT Multi-unit Abutment WP ex-hex (Brånemark)

² Compatible ONLY with Multi-unit Abutments WP ex-hex (Brånemark)

Universal Base non-engaging: clinical screw included.

On1™ Universal Base

Single unit restoration – Engaging

35
Uni



On1 Universal Abutment – Engaging

Connection	On1 Base		
Platform	NP	RP	WP
Article nr. collar 0.3 mm	38707	38709	38711
Article nr. collar 1.25 mm	38708	38710	38712
Post height	3.98 mm	3.98 mm	3.98 mm
Diameter	Ø 4.775 mm	Ø 4.775 mm	Ø 4.775 mm
Screw	38750	38750	38750
Screwdriver	Unigrip	Unigrip	Unigrip
Lab Screw	38751	38751	38751
IOS scanbody & Healing cap (4.5 mm / 5.0 mm / 6.0 mm)	38735 / - / 38734	38737 / - / 38736	38739 / 38738 / -
IOS / Desktop analog	38740	38741	38742
Insertion tool	38201	38203	38204
Desktop scanbody DTX	300437	300437	300437
Plaster model implant analog	38731	38732	38733

Delivered with the corresponding On1 Prosthetic screw and On1 Burn-out coping.

On1™ Universal Base

Multiple unit restoration – Non-Engaging

35
Uni



On1 Universal Abutment – Non-engaging

Connection	On1 Base		
Platform	NP	RP	WP
Article nr. collar 0.3 mm	300361	300363	300365
Article nr. collar 1.25 mm	300362	300364	300366
Post height	3.98 mm	3.98 mm	3.98 mm
Diameter	Ø 4.775 mm	Ø 5.3 mm	Ø 6.515 mm
Screw	38750	38750	38750
Screwdriver	Unigrip	Unigrip	Unigrip
Lab Screw	38751	38751	38751
IOS scanbody & Healing cap (4.5 mm / 5.0 mm / 6.0 mm)	38735 / - / 38734	38737 / - / 38736	38739 / 38738 / -
IOS / Desktop analog	38740	38741	38742
Insertion tool	38201	38203	38204
Desktop scanbody DTX	300437	300437	300437
Plaster model implant analog	38731	38732	38733

Delivered with the corresponding On1 Prosthetic screw and On1 Burn-out coping.

Screwdrivers

OmniGrip™

OmniGrip™ Screwdriver Machine

20 mm	37379
25 mm	37380
30 mm	37381
35 mm	37382



OmniGrip™ Screwdriver Manual

20 mm	37376
28 mm	37377
36 mm	37378



OmniGrip™ Mini

OmniGrip™ Mini Screwdriver Machine

20 mm	300852
28 mm	300853
36 mm	300854



OmniGrip™ Mini Screwdriver Manual

20 mm	300855
36 mm	300856



UniGrip™

UniGrip™ Screwdriver Machine

20 mm	29151
25 mm	29152
30 mm	29153
35 mm	29154



UniGrip™ Screwdriver Manual

20 mm	29148
28 mm	29149
36 mm	29150



Tools

Elos Accurate™ Analog insertion tool and IO driver

Elos Accurate™ Analog Insertion Pin	PMA-AIP-2
Elos Accurate™ Analog Insertion Screw	PMA-AIS-5
Elos Accurate™ Analog Pliers	PMA-AIT-1
Elos Accurate™ IO driver – mini	C18512
Elos Accurate™ IO driver – short	C13485
Elos Accurate™ IO driver – long	C13048

Compatible softwares



Other tools – Manual Torque Wrench

Manual Torque Wrench Prosthetic	29165
Manual Torque Wrench Adapter	29167



Manual Torque Wrench Prosthetic Nobel Biocare N1™	301082
---	--------



For availability of compatible software and download of local production DME's, please go to <https://www.nobelbiocare.com/local-production>