


# SLIM FLUSH CUT & HOLD DISTAL END CUTTER

| *O-Ring Design, 678-114*



Hu-Friedy's inserts ensure the cutter flushly cuts wire distally to the buccal tube or bracket. The blue silicone ring holds the cut wire until released. The slim design of the cutter increases intraoral visibility, reduces debondings and is more comfortable for the patient.

The 678-114 is designed for all wire types, up to a maximum of .021" x .025" (.53mm x .64mm).

## POINTS OF PERFORMANCE

### **Crafted with Immunity Steel®**

A proprietary blend of Surgical Grade Stainless Steel for a long-lasting cutting edge plus corrosion resistance.

### **Cuts Flush**

Distal wire is cut flush to brackets for optimal patient comfort.

### **Slim Design**

Slender working end improves intraoral visibility, reducing debondings.

### **Blue Ring**

Brightly colored medical grade silicon ring cushions and holds the cut wire.

10 blue silicon rings included in package. Replacement rings can be purchased.

| **678-114R50** 50 PACK OF RINGS or | **678-114R150** 150 PACK OF RINGS

**Get a demonstration on how to place the blue ring.**

Search YouTube for "Hu-Friedy 678-114 ring placement."

| **678-114**

SLIM FLUSH CUT & HOLD DISTAL END CUTTER

Visit us online at [Hu-Friedy.eu](http://Hu-Friedy.eu)

©2018 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-085GB/0818

How the best perform

