


# HU-FRIEDY HD BLACK LINE MIRRORS

| Proven to Reduce Glare up to 80%\*



See your patients - not your instruments. Glare interferes with visibility. Reducing this disruption to your visibility allows you to see more detail and offer more accurate patient care. HD Black Line Mirrors are engineered to optimize clinical outcomes by delivering superior visibility throughout any dental procedure. The reduced glare properties of the Diamond Like Carbon (DLC) coating and clarity of the HD mirror glass facilitates quicker and more accurate visibility intraorally.



## REDUCES GLARE UP TO 80%\*

“ I wear loupes with a LED light and noticed a significant difference in the amount of glare reflecting back into my field of vision during patient care, allowing me to see more detail. I realized, at the end of my clinical day, I had less eye strain when I used the HD Black Line Mirrors. ”

– Tami Wanless, RDH, MED

### DIAMOND LIKE CARBON (DLC) COATING

The **durable black matte coating** reduces the glare compared to a standard metal mirror head and handle. This helps to **reduce strain and fatigue** as the user does not need to adapt their viewing position due to unwanted shine.

### ENHANCED CONTRAST & VISUAL ACUITY

The **black color** of the coating creates a distinct contrast between the instrument, the tooth and/or the surrounding tissue allowing for easy identification intraorally.\*

### ERGONOMIC HANDLE

The large diameter handle increases control, provides maximum comfort and reduces hand fatigue, creating an ergonomically friendly handle option.

### SUPERIOR BRILLIANCE & COLOR

HD Mirror provides **superior brilliance and color** for quicker and more **accurate visibility** of the mouth.†

- 113%<sup>†</sup> reflection factor for exceptional image clarity
- 38,5%<sup>†</sup> brighter than rhodium coated mirror glass
- 50%<sup>†</sup> brighter than other front surface mirror glass

Visit us online at [Hu-Friedy.eu](http://Hu-Friedy.eu)

©2019 Hu-Friedy Mfg. Co., LLC. All rights reserved. \*When comparing the Hu-Friedy DLC coated mirror head and handle to the Hu-Friedy non-coated stainless steel mirror head and handle.  
† Data on file and available upon request.

How the best perform




# ENHANCE YOUR SET-UPS WITH HD BLACK LINE MIRRORS


DIAGNOSTIC




**Hu-Friedy's new HD Black Line Mirrors are not just for Hygiene Exam procedures. Complete your Kit no matter the setup!\***




General Exam




Basic Hygiene



Exam Perio



Restorative



Endodontic



Surgical

## PARTCODES & DESCRIPTIONS

### CONE SOCKET THREAD

### EUROPEAN THREAD



#### MIR5DSHDB

#5 Double-sided HD Mirror, Cone Socket, Black Line, Single Pack



#### MIR5HDB

#5 Single-Sided HD Mirror, Cone Socket, Black Line, Single Pack



#### MIR5DSHDBE

#5 Double-sided HD Mirror, European, Black Line, Single Pack



#### MIR5HDBE

#5 Single-Sided HD Mirror, European, Black Line, Single Pack

#### MH6B

Cone Socket Mirror Handle, Satin Steel, Black Line



#### MHE6B

European Mirror Handle, Satin Steel, Black Line

Manufacturer: Hu-Friedy Mfg. Co., LLC | 3232 N. Rockwell Street | Chicago, IL 60618 | USA  
Customer Care Department: Hu-Friedy Mfg. Co., LLC.  
Free Call +800 HUFRIEDY (00800 48 37 43 39) | Free Fax 00800 48 37 43 40  
E-Mail: info@hu-friedy.eu | Website: Hu-Friedy.eu

©2019 Hu-Friedy Mfg. Co., LLC. All rights reserved. \*HD Black Line Mirrors sold separately. HFL-285GB/0919

How the best perform

