

KATANA™ Zirconia



New 90 Minutes Sintering
STML Shade Availability Expansion

UTML *Ultra Translucent Multi Layered*

STML *Super Translucent Multi Layered*

ML *Multi Layered* **HT*** *High Translucent*



The logo "MULTI LAYERED" expresses original production technology of KATANA™ multi-layered materials.

TOTAL SINTERING TIME: 90 MINUTES

KATANA™ Zirconia, the innovative multi-layered zirconia, known for specialized fabrication of Full Contour Zirconia prosthetics is now speed capable. Fast sintering furnaces can now be used with KATANA™ Zirconia in the speed sintering process*¹. You can successfully finish sintering KATANA™ Zirconia in approximately 90 minutes*².


*¹ Only for single unit *² The total sintering time may be affected depending on: furnace condition; surrounding temperature; voltage; and electrical current, etc.

Fast sintering furnaces



AUSTRUMAT 674i (DEKEMA Dental-Keramiköfen GmbH)
inFire HTC speed (Dentsply Sirona)
Sintra Plus (Shenpaz Dental Ltd.)

Sintering schedule



EXPANDING THE SHADE AVAILABILITY OF STML

STML exhibits an excellent balance of translucency and desired mechanical properties. Now available in 9 expanded shades for completion of a wider range of esthetic restorations.



Printed color can be slightly different from actual color.

YOUR CONTACT

Kuraray Europe GmbH
BU Medical Products
Philipp-Reis-Straße 4
65795 Hattersheim am Main

Phone DACH
Phone Rest of Europe
Fax
E-Mail
Website

+49 (0) 69-305 35 835
+49 (0) 69-305 85 980
+49 (0) 69-305 98 35 835
dental.eu@kuraray.com
www.kuraraynoritake.eu



Kuraray Noritake Dental Inc.

300 Higashiyama, Miyoshi-cho, Miyoshi, Aichi 470-0293, Japan

*KATANA™ is a trademark of NORITAKE CO., LIMITED.