



PRIRODNA LJEPOTA ZA VAŠE PROTETSKE RADOVE

Saznajte više o KATANA™ ZIRCONIA porodici



KATANA™ Zirconia HTML High Translucent Multi Layered

Daje čvrste i estetske restauracije.
Struktura od cirkona poput dentina? Da!

S KATANA Zirconia HTML svjetlo reflektira prirodnu opaknost. Nanošenje keramike na površinu strukture daje najprirodniji rezultat. Ovaj cirkon dostupan je u paleti prirodnih boja. Korištenje ovih diskova idealno je za rad bezmetalnih struktura. HTML diskovi pogodni su za rad krunica i mostova u punoj formi.



KATANA™ Zirconia STML Super Translucent Multi Layered

Idealan odnos između translucencije i jačine.
Da li mogu cirkonski protetski radovi u punoj formi izgledati prirodno? Da!

S KATANA Zirconia STML blokovima svjetlo se prenosi u incizalno područje i zaustavlja u cervikalnom području. Na ovaj način imate fleksibilnost izrade lijepih protetskih radova u punoj formi, bez obzira na boju nadogradnje. Ovaj proizvod je prvi na tržištu sa karakteristikom cirkona s bojom i prirodnom translucencijom i gradacijom boje. Rezultat rada s ovim diskovima su izrazito estetski i čvrsti protetski radovi.



KATANA™ Zirconia UTML Ultra Translucent Multi Layered

Translucencija slična caklini
Radite li ljuskice prirodnog izgleda i ostale protetske radove u prednjoj zoni od cirkona? Da!

KATANA Zirconia UTML ima izrazitu transparentnost, zahvaljujući čemu se postiže prirodna boja zuba te rezultat prirodnosti koju želimo postići. Ovaj cirkon ima karakteristike najveće transparentnosti na tržištu uz paletu raspoloživih prirodnih boja zuba. Proizvod je izvrstan i jedinstven!

MULTI
LAYERED

KATANA Zirconia
Multi-Layered serija ima
tri važne karakteristike:
ekonomičnost, brzina i
preciznost. Digitalna tehnologija
daje rješenje i u estetskoj zoni.

14 mm

1.260 Kn
+ PDV 5%

18 mm

1.590 Kn
+ PDV 5%

22 mm

1.790 Kn
+ PDV 5%



BOJE

HTML	NW	A1	A2	A3	A3.5	A4	B1	B2	B3
	C1	C2	C3	D2	D3				

STML	NW	A1	A2	A3	A3.5	A4	B1	B2	B3
	C1	C2	C3	D2	D3				

UTML	A1	A2	A3	A3.5	A4	B1	B2	B3	B4
	C1	C2	C3	C4	D2	D3	D4		
	ENV	EA1	EA2	EA3					

30-75% jači od e.maxa
Test savijanja u 3 točke



CZR™ FC PASTE STAIN

Inovativni materijal za jednostavnu karakterizaciju restauracija od punog cirkona.

Dostupno u 27 boja, velika mogućnost reprodukcije boja. Široka paleta plavih i sivih boja daje izvrsnu mogućnost transparentnog efekta u incizalnom području. Dva oblika glazure (Glaze i Clear Glaze) olakšavaju mogućnost transparentnosti kroz cijelu krunicu.



Enamel	Dentin	Cervical	Incisal Effect	Body Effect	Pink Tissue
Glaze	A+	Cervical 1	Light Gray	Mamelon Orange 1	Salmon Pink
Clear Glaze	B+	Cervical 2	Dark Gray	Mamelon Orange 2	Pink
Value	C+	Cervical 3	Grayish Blue	Earth Brown	Red
Fluor	D+		Blue	Reddish Brown	
			Black	Yellow	
			White	Green 1	
				Green 2	



CZR™ FC Paste Stain Full Kit

Sadrži: Glaze, Clear Glaze (5g cad.); Light Gray, Dark Gray, Grayish Blue, Blue, Cervical1, Cervical2, Cervical3, A+, B+, C+, D+, White, Mamelon Orange1, Mamelon Orange2, Earth Brown, Reddish Brown, Value, Black, Green1, Green2, Yellow, Salmon Pink, Pink, Red, Fluor (3g cad.); ključ boja; ES Liquid (10ml)

CZR™ FC Paste Stain Full Kit

(cod. 104-5590-EU)

7.950 Kn
+ PDV 5%



CZR™ FC Paste Stain Basic Kit

Sadrži: Clear Glaze (3g.); Dark Grey, Grayish Blue, Cervical2, Earth Brown, A+, White, Salmon Pink (3g.); ključ boja; ES Liquid (10ml)

CZR™ FC Paste Stain Basic Kit

(cod. 104-5870-EU)

2.490 Kn
+ PDV 5%

KATANA™ ZIRCONIA BLOCK STML


Multi-layered KATANA STML Zirconia Block sa 4 razine cirkona u gradiranim bojama.

Proizvod namijenjen za rad u ordinaciji u prirodnim bojama – efikasniji i brži način rada. Jačina 763 MPa!
15 minuta obrada + 18 minuta sintranje

Dostupno u 3 veličine: 12Z, 14Z i 14Z L
Dostupno u bojama NW, A1, A2, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, D2 i D3.
Dostupno za CEREC.

Indikacije:

- Krunice u prednjoj i stražnjoj zoni
- Ljuskice
- Inlej - onlej
- Most 3 članova



		a (height)	b (width)	c (length)
12Z	Before sintering	15.3 mm	19.2 mm	20.2 mm
	After sintering	12.2 mm	15.4 mm	16.2 mm
14Z	Before sintering	17.8 mm	19.2 mm	20.2 mm
	After sintering	14.2 mm	15.4 mm	16.2 mm
14Z L	Before sintering	17.8 mm	19.2 mm	40.0 mm
	After sintering	14.2 mm	15.4 mm	32.0 mm

12Z (5 kom)
990 Kn
+ PDV 5%

14Z (5 kom)
1.050 Kn
+ PDV 5%

14Z L (5 kom)
1.150 Kn
+ PDV 5%



KATANA™ AVENCIA™

Kompozitni blokovi
Hibridna keramika





Dostupno u 2 veličine:
12 (jedna krunica) i 14L (most)
Boje: A2, A3 i A3.5
Dostupno za: CEREC i UNIVERSAL nastavke

Indikacije:

- Inlej - onlej
- Premolarne krunice
- Krunice u prednjoj i stražnjoj zoni
- Most 3 članova

	a (height)	b (width)	c (length)
12	10 mm	12 mm	15 mm
14L	14.5 mm	14.5 mm	18 mm



12 Universal
(5 kom)

790 Kn
+ PDV 5%

14L Universal
(5 kom)

840 Kn
+ PDV 5%

12 Cerec
(5 kom)

870 Kn
+ PDV 5%

14L Cerec
(5 kom)

920 Kn
+ PDV 5%

PANAVIA Cementi



PANAVIA V5

Panaviu V5 odlikuje iznimna čvrstoća i naglašena estetska komponenta!

Panavia V5 je najjači cement ikad razvijen, a dostupan je u 5 nijansi. Koristi se za cementiranje krunica, mostova, ljuskica izrazito je pouzdan za nagoradnje na implantatima. Postupak započinjete s tretiranjem zuba s PANAVIATMV5 primerom, nastavljate s tretiranjem protetskog rada s CLEARFILTM CERAMIC PRIMER PLUS, potom cementirate - i to je to. PANAVIATMV5AUTOMIX CEMENT olakšava rukovanje i apliciranje uz jednostavno uklanjanje viška cementa. Automix šprica omogućava brz, jednostavn i precizan rad - a konzistentno optimiziran cement omogućuje da postavite protetski rad na način koji želite.

PANAVIA™ V5 Professional Kit

Tooth Primer (2 ml), Ceramic Primer Plus (2 ml), Panavia V5 Pasta u 5 boja (2.4 ml/4.2g), probne paste u 5 boja, ukupno 12 ml u 5 boja.

2.790 Kn
+ PDV 5%

PANAVIA™ V5 pasta

pasta 4,6 ml + 20 mješalica

890 Kn
+ PDV 5%

PANAVIA™ V5 Uvodni set

Boja Universal (A2) ili Clear Tooth Primer (2 ml), Ceramic Primer Plus (2 ml), Panavia V5 Pasta (2.4 ml/4.2g)

950 Kn
+ PDV 5%



PANAVIA SA CEMENT UNIVERSAL

Samoadhezivni kompozitni cement

Jedinstveni i inovativni kompozitni cement, samoadhezivni, indiciran za sve vrste materijala (keramiku, litij disilikat, kompozit) bez potrebe silaniziranja. Uz jedinstveni MDP monomer u formuli je dodana komponenta LCS (LONG CHAIN SILANE) koja osigurava snažnu kemijsku vezu na sve vrste površina. To je samojetkajući cement koji ne treba korištenje primera i može se koristiti u svim slučajevima - suhom i vlažnom okruženju; kod cementiranja zirkona, kompozita, litij disilikata (E-max krunica). Zahvaljujući MDP monomeru cement ima funkciju samoadhezije na bazi kemijske veze. U cement je uključena jedna nova komponenta LCSi MONOMER (LONG CHAIN SILANE MONOMER). To je jedinstvena formula koju nema niti jedan drugi cement - što znači da nije potreban silan za krunicu (kao Ceramic Primer ili slično).

PANAVIA™ SA Cement Universal

Boja A1 ili translucent, 1 x tubica 8,2 g, sa priborom

780 Kn
+ PDV 5%

VALUE
PACK
2+1

PANAVIA™ SA Cement (3kom)

Boja A1 ili translucent, 1x tubica 8,2 g, sa priborom

1.560 Kn
+ PDV 5%

KATANA™ TWIST DIA

3 koraka za blistav sjaj zirkona



KATANA™ TWIST DIA

3 kom

255 Kn
+ PDV 25%

KATANA CLEANER

Prvi univerzalni intra i ekstra oralni restaurativni cleaner površine.

KATANA™ Cleaner je prvo univerzalno sredstvo za čišćenje svih površina koji ima pH 4.5 što omogućuje ne samo ekstraoralno već i intraoralno čišćenje.

Indikacije:

- Intraoralna primjena: zub, bataljak, korijenski kanal
- Ekstraoralna primjena: cirkonske strukture, litij disilikat, keramika, materijali od kompozita CAD-CAM krunice, kompozit, metali (legure plemenitih / neplemenitih metala), kolčiči (kompozitni i metalni)



DEKONTAMINACIJA ZA 10 SEKUNDI

KATANA™ Cleaner

4 ml

198 Kn
+ PDV 25%



Clearfil Ceramic Primer plus

Novi Clearfil Ceramic Primer se odlično slaže sa svim vrstama keramike. Osim već dokazanog monomera MDP, koji osigurava snažnu vezu, jednokomponentni primer također sadrži etanol i silan agent MPS, koji osigurava snažnu vezu na sve materijale od kojih se izrađuju estetske restoracije. Kroz optimalnu kombinaciju ta tri elementa, Clearfil Ceramic Primer održava izvrsnu dugotrajnu adheziju.

Clearfil Ceramic Primer plus

4 ml

690 Kn
+ PDV 25%

EX-3™ Keramika za metal

*10 x snažniji na kobalt-kromu
Adheziv koji poboljšava vezu između keramike i plemenitih legura, na bazi krom-kobalta



EX-3™ N Color Intro Kit

Sadrži: Paste Opaque NA2 (6g), Paste Opaque NA3 (6g), Body NA2B (10g), Body NA3B (10g), Enamel E2 (10g), Enamel E3 (10g), Luster LT1 (10g), Luster Creamy Enamel (10g), Luster CCV-1 (10g), Translucent ADT (10g), Translucent TX (10g), Forming Liquid (10ml)

EX-3™ N Color
Intro Kit
(cod. 102-0050-EU)
1.429 Kn
+ PDV 5%

EX-3™ NP Bonder
6 g
(cod. 102-4532-EU)
241 Kn
+ PDV 5%

EX-3™**
10 g
94 Kn
+ PDV 5%
* opcije boja
Opacius Body OBA3
Clear Cervical CCV-2
Speed Enamel S1
Speed Enamel S2
Speed Enamel S3
Speed Enamel S4

CZR™ Keramika za cirkon



CZR™ Luster CCV Kit

Sadrži: Lt0 (10G), Lt1 (10G), T Blue (10G), Aqua Blue 1 (10G), Aqua Blue 2 (10G), Lt Natural (10G), Lt Super Gray (10G), Creamy Enamel (10G), Sun Bright (10G), Incisal Aureola (10G), Creamy White (10G), Lt Yellow (10G), Ccv-1 (10G), Ccv-2 (10G), Ccv-3 (10G), Ccv-4 (10G), Ccv-2 (10G), Ccv-3 (10G), Ccv-4 (10G).

CZR™ Luster CCV Kit
(cod. 104-1790-EU)
2.322 Kn
+ PDV 5%



CZR™ Internal Stain Complete Kit

Sadrži: White (3G), Incisal Blue 1 (3G), Incisal Blue 2 (3G), Mamelon Orange 1 (3G), Mamelon Orange 2 (3G), Cv-1 (3G), Cv-2 (3G), Cv-3 (3G), Earth Brown (3G), Reddish Brown (3G), Fluoro (3G), Salmon Pink (3G), Red (3G), A+ (3G), B+ (3G), C+ (3G), D+ (3G), Gray (3G), Bright (Dilution) (3G), Is Liquid (10 ml), Ključ Boja.

CZR™ Internal Stain Complete Kit
(cod. 104-1840-EU)
3.600 Kn
+ PDV 5%



CZR™ Trial Kit SS A2

Sadrži: OBA2 (5g), A2B (5g), E2 (5g), LT1 (5g), IS Liquid (10ml), Forming Liquid (10ml), CZR PRESS SBS SSA2 (3g).

CZR™ Trial Kit SS A2
(cod. 104-5530-EU)
300 Kn
+ PDV 5%



PERFECTA 900

Laboratorijski motor za cirkon
Dolazi u opciji sa ugrađenom turbinom,
a dodatno se može spojiti i motor LA-66
i time dobiti kompletno rješenje za lab

- motor 7.8 Ncm / 1.000 - 50.000 rpm
- snažna turbina 0,7 Ncm / 5.000 - 100.000
- 2 verzije hlađenja: integrirani tank i vanjski spoj
- 2 verzije uređaja: stolna i koljenska
- Air-jet funkcija
- 24 mjeseca garancije

PERFECTA 900
16.400 Kn
+ PDV 25%

Motor LA-66
za ugradnju na
PERFECTA 900
8.025 Kn
+ PDV 25%



PERFECTA 600

Laboratorijski motor za cirkon
sa motorom LA-66

- motor 7.8 Ncm / 1.000 - 50.000 rpm
- Air-jet funkcija
- 3 verzije: stolna, koljenska i nožna
- 24 mjeseca garancije

PERFECTA 600
15.200 Kn
+ PDV 25%



novodent

Uvoznik & distributer

Novodent d.o.o. • info@novodent.hr • www.novodent.hr

MATULJI Cesta dalmatinskih brigada 7 • 51211 Matulji • Tel: 051 701 645

ZAGREB Zavrtnica 3 • 10000 Zagreb • Tel: 01 457 7846

SPLIT Tel: 091 2 704 712